

Optical Imaging meets MS Imaging: Integrating a Dual-Beam Scanning Electron Microscope and a SIMS into one Device for Simultaneous Imaging



University of Wuppertal, Germany
Institute for Pure and Applied Mass Spectrometry

Alexander Laue¹, Albrecht Glasmachers¹, Albrecht Brockhaus¹, Michel Aliman²


¹University of Wuppertal, Wuppertal, Germany, ²Carl Zeiss NTS GmbH, Oberkochen, Germany

Introduction

Overview:

Optical imaging using a Scanning Electron Microscope (SEM) is a well-known tool in the nanotechnology, semiconductor industry and material science which allows for nanometer resolution.


For advanced imaging capabilities, sample preparation and surface processing a Focused Ion Beam (FIB) is used. The FIB produces a small amount of secondary ions which can be used for high laterally and spatially resolved mass spectrometry imaging.



Transfer efficiency of the SIMS Device

Primary ion beam current: 100 pA
Secondary ion beam current: approx. 200 fA
⇒ 125.000 ions in 100 ms


Electrostatic extraction lens [1]:
• broad kinetic energy distribution (Thompson)



IonGuide [2,4]: Linear Quadrupole
• axial ion transfer (segmented dc-gradient)
• ion cooling with helium as damping gas
• ion bunching
• mass filter

WaveGuide[3,4]: Miniature Quadrupole
• accumulating secondary ions
• sequential ion transfer
• pressure stage

Electrical ion resonance cell [5,6]:
• FFT Mass analyzer



Extraction efficiency: 12 %
(simulation & estimation)
approx. 15.000 ions

transfer efficiency: (IonGuide): > 40 %
approx. 6.000 ions

transfer efficiency: (WaveGuide): > 70 %
approx. 4.200 ions
kinetic energy of ion: < 0,5 eV FWHM

Overall transfer efficiency inside the SIMS: 28 %

Conclusions

- The combination of SIMS, FIB and SEM opens doors to new applications and methods which are impossible with separate devices
- A full functional SIMS for Scanning Electron Microscopes has been built and evaluated
- The SIMS has been tested with different anorganic and structured samples
- High transfer and cooling efficiencies of sputtered secondary ions have been shown
- High sensitivity of the mass analysis has been shown at FIB currents down to 10 pA
- To avoid surface defects low FIB currents and low sputter rates must be used
- The FFT mass analyzer allows for multiple measurements of the same ion population with higher sensitivity
- Scanning speeds up to 10 Hz have been demonstrated

Applications & Requirements

Mass spectrometry controlled material processing:

- Material dependent sputtering (e.g. end-point detection) - only remove those parts of a target which consist of a specific material
- Advanced layer thickness analysis - detect different layers whose material contrasts are too low for optical detection

Requirements:



- Real-time mass spectra generation
- Short measurement time per spectrum

Optical controlled mass spectrometry:


- Optical feedback of sputtered spot (simultaneous imaging)
- Find and identify spots of interest without surface defect (SEM) and create selective MS spectra with minimal surface defect (FIB+SIMS)
- Analyze optically identified spots of interested (fractional mapping)

Requirements for MS imaging:


- Fast sampling rate per MS-Pixel: 100 ms
- High lateral and axial resolution
- Low target defect ⇒ high sensitivity




Measurement Results




Ion transient:
gallium and molybdenum ion signal




Spectrum of Molybdenum:
300 pA FIB current (Gallium) @ 30 keV,
100 ms accumulation time




High sensitivity:
10 pA FIB current
300 ms accumulation time



Sample target:
TEM grid (Cu) filled with
conducting silver



Resulting spectrum of TEM grid:
700 pA FIB current (Gallium) @ 30 keV,
200 ms accumulation time



Target defect:
Sputtering of a few seconds leads to massive implantation of gallium

References

- [1] Prellszas, D.; Aliman M.; Mantz, H.; Laue, A.; Brockhaus, A.; Glasmachers, A.; Optimization of the collection efficiency of secondary ions for spatially resolved SIMS in Crossbeam devices 12th International Seminar on Recent Trends in Charged Particle Optics and Surface Physics Instrumentation, Bmo (Czech Republic), 2010
- [2] Laue, A.; Glasmachers, A.; Brockhaus, A.; Aliman, M.; Highly efficient and compact Quadrupole Ion Transfer Guides in planar technology for space-limited applications 18th International Conference on Secondary Ion Mass Spectrometry, Riva del Garda (Italy), 2011
- [3] Glasmachers, A.; Laue, A.; Brockhaus, A.; Aliman, A.; Transfer Efficiency and Timing Performance Measurements of Multipole Ion Guides and Ion Wave Guides Constructed with Planar Technologies 59th ASMS Conference, Denver (USA), 2011
- [4] Glasmachers, A.; Laue, A.; Brockhaus, A.; Püwey, A.; Aliman, M.; Planar technologies for optimized realizations of quadrupole ion guides and quadrupole ion wave guides 58th ASMS Conference, Salt Lake City (USA), 2010
- [5] Laue, A.; Glasmachers, A.; Brockhaus, A.; Aliman, A.; Mantz, H.; Metrological Characterization of a Sensitive Secondary Ion Mass Spectrometer for Electron Microscopes to Combine Optical/Structural and Analytical Imaging 59th ASMS Conference, Denver (USA), 2011
- [6] Laue, A.; Glasmachers, A.; New Design of a Compact Fourier-Transform Quadrupole Ion Trap for High Sensitivity Applications 57th ASMS Conference, Philadelphia (USA), 2009
- [7] Patent application EP11152379.1 -1232: Apparatus for focusing and for storage of ions and for separation of pressure areas, 2011
- [8] Patent application EP11152420.3-2208: Apparatus for transmission of energy and/or for transportation of an ion as well as particle beam device having an apparatus such as this, 2011